



Product Portfolio

Science. Results. Trust.

SYNERON  CANDELA®

INNOVATION, BROUGHT TO YOU BY SYNERON CANDELA.

For decades, Syneron Candela has been a trusted name for dermatologists, plastic surgeons and aesthetics specialists the world over. We have engineered the systems that serve as the basis for building profitable practices—reliable solutions that use proven energy-based technologies to treat a vast array of conditions, as well as the full spectrum of skin types.

AND THAT'S JUST THE BEGINNING.

Syneron Candela is dedicated to developing technology that generates even more impressive results and extends your power to treat even more conditions. We pioneered the use of pulsed dye lasers to treat a range of conditions, safely and effectively. Under the

leadership of our founder Shimon Eckhouse, we developed patented elos® (electro-optical synergy) technology, combining bi-polar radio frequency with light-based energy to offer safer, more comfortable treatments that deliver incredible results.



2000: SYNERON

Develops elos technology, led by Shimon Eckhouse, Inventor of Intense Pulsed Light (IPL) technology



1970: CANDELA



Develops the first laser based on the theory of selective photothermolysis

1985

Candela begins testing on the first Pulsed Dye Laser, now the industry's gold standard for safety and efficacy

1993

The first Q-switched nanosecond technology for the treatment of unwanted tattoos and pigmented lesions is launched in the market

1998

Candela's Alexandrite laser is cleared by the FDA for permanent hair reduction. It's now the most-used hair removal system worldwide

2002

Candela launches Smoothbeam 1450 nm diode for the treatment of acne vulgaris, acne scars and periorbital wrinkles



893 patents granted and pending



29% global market share¹



The choice of **30,000** physicians worldwide




700

peer-reviewed articles



100

clinical studies in progress



“ We are committed to industry leadership by focusing on continuous technological innovation that helps improve and change lives.”

Shimon Eckhouse

2004

Syneron is first publicly traded. The first system incorporating RF and IR for the temporary improvement of cellulite

2007

Syneron launches first light-based device for the treatment of cellulite

2008

Syneron launches first RF technology for aesthetic applications

2009

Syneron begins partnership with Procter & Gamble

2010
Syneron acquires Candela

2012
Syneron acquires UltraShape and introduces the first pulsed ultrasound for fat reduction

2014 – 2016
Syneron Candela introduces PicoWay picosecond laser, with fractional capabilities and three wavelengths

The Syneron Candela Portfolio

Aesthetic solutions for both the face and body



PicoWay® - The Complete PicoSecond Platform

3 wavelengths – 532nm, 785nm & 1064nm wavelengths in one laser system to treat a wide variety of pigmented lesions and all tattoo colors & types, including recalcitrant tattoos.

NEW Resolve™ dual wavelength picosecond fractional module for treating pigmentation, skin irregularities and signs of aging.

High peak power – For effective energy over the broadest range of spot sizes.

Shortest picosecond pulses – 40% shorter picosecond pulses mean effective, yet lower energies for less risk of side effects.



Gentle Pro Series™ - Integrated aesthetic treatment workstation

Fast and powerful, dual-wavelength, multi-spot laser platforms deliver treatment efficacy, flexibility and safety.

Treat unwanted hair, pigmented lesions, vascular lesions and other indications comfortably and effectively on all skin types.

Choice of cooling systems: Chilled air or cryogen based DCD™ for convenient and consistent epidermal protection.

Available as cost-effective, upgradable single wavelength 1064 nm and 755 nm systems.

20, 22 & 24 mm treatment spots are the largest in the market, enabling shorter treatment times and improved patient comfort.



CO₂RE® - The Complete CO₂ Solution

Offers proven solutions for skin resurfacing, wrinkle and scar reduction, pigmentation, minor dermal incisions and excisions, and now vaginal treatments.

A powerful CO₂ in a compact device - 6 treatment modes for optimal tissue impact and total flexibility.

Unique Fusion mode treats both superficial and deep skin layers at the same time.

In-board laser scanners mean that CO₂RE's high quality handpieces are ergonomically shaped and incredibly light.



Partner 360 is part of Syneron Candela's commitment to providing easy access to the marketing tools you need to grow your practice. This user-friendly customer portal and eStore enables you to quickly and easily download customizable practice marketing materials and a wide range of marketing resources 24 hours a day, 7 days a week.



eIos Plus® - Complete aesthetic work station

The ultimate multi-application platform, with 15 different applicators for the most in-demand treatments.

eIos Plus uses eIos® technology (a combination of optical energy with bi-polar radio frequency) to enable safe and effective treatments on all skin types.

NEW trinit Plus - 3 Treatments in 1 Session - treat color correction, contouring, wrinkles and textures.

NEW Motif Vantage, Fast Diode Hair Removal, has the largest spot size available for fast treatments.

Intelligent Feedback System (IFS) ensures precise delivery of RF energy regardless of variations in skin impedance for consistent results and added safety.

Intuitive user interface with easy-to-use guided treatment modes.



UltraShape Power - Body shaping via fat cell destruction

The new **UltraShape Power** comes with **20% increase** in acoustic energy and **32% average fat layer reduction**.

UltraShape focused ultrasound technology has been clinically proven to be safe and effective with multiple peer-reviewed papers and over **400,000 treatments** performed worldwide.

Coverage

Over 50% more body area treated per hour than competitors.

Customization

Tailored to the contours of the patient's body.

Cost

Higher profitability than any other non-invasive fat-reducing device.

Combination Therapy

UltraShape VelaShape for fat destruction and tightening & cellulite.

Comfort

Mechanical effect, not thermal, for complete patient comfort with no bruising, swelling or downtime.



VelaShape III with Guided Mode - The art of body contouring

Top-performing **Body Contouring Solution**, Now Better Than Ever.

Improve outcomes for **body contouring, circumference reduction, cellulite improvement, and skin tightening**.

Excellent results in just **1-3 sessions**.

Effective complementary therapy-

Boosts results from other fat destruction procedures, such as UltraShape and liposuction.

Now, **VelaShape III with Guided Mode** offers an enhanced treatment experience – for both patients and operators with simplified operation.

With more than **1 million VelaShape procedures** performed yearly, VelaShape III has proven time and its ability to deliver significant results.

A complete portfolio of innovative technologies for aesthetic therapies



Vbeam® - Best-in-class pulsed dye laser

The gold standard in vascular treatments.

More published studies than any other aesthetic laser on the market.

Specifically cleared for all ages - treat infants, children and adult patients.

Safe and clinically proven - an excellent addition to any aesthetic practice.

Expanding aesthetic applications - reduces bruising associated with injectables and other cosmetic procedures.



Alex TriVantage® - Pigmented lesion & tattoo removal treatments

The gold standard treatment for benign dermal and epidermal pigmented lesions.

Unique patented laser-pumped laser handpieces offer 755, 532 and 1064 nm treatment capabilities.

Treats a wide range of tattoos - no other platform treats such a wide variety of tattoo colors and other pigmented lesions at variable depths.

The perfect complement to any aesthetic practice.



Profound® - 1st device to create all 3 skin fundamentals

Creates elastin, collagen and hyaluronic acid, profoundly.

Creates dermal volume, naturally.

Yields 100% response¹ for total patient satisfaction.

Treats jawline, jowls, wrinkles and other textural considerations.

NEW SubQ Module to treat adipose tissue and septae.

Yields predictable results with complete, real time control over depth, temperature and time of treatment.

Over 10 million treatments performed annually with our products

¹Based on results of a clinical study, measuring improvement in Fitzpatrick Wrinkle Scale at 3 months compared to baseline based on independent review of photographs.



Partner 360 is part of Syneron Candela's commitment to providing easy access to the marketing tools you need to grow your practice. This user-friendly customer portal and eStore enables you to quickly and easily download customizable practice marketing materials and a wide range of marketing resources 24 hours a day, 7 days a week.



eTwo® - Complete facial rejuvenation

The original fractionated RF platform lightweight, compact, applicator table-top device.

Sublime uses elōs technology to remodel the deep dermal layer of the skin by stimulating new collagen growth.

Sublative increases collagen by placing fractionated heat energy effectively into the upper dermis with minimal epidermal disruption.

Quick treatments - typically 20 to 40 minutes depending on area treated.



eLase® - Fast, pain free, effective hair removal

Fastest diode laser platform- delivers safe and effective treatments on the widest variety of skin and hair types.

Motif LHR® Hair Removal includes the high frequency and lower energy Motif mode of operation for fast and virtually pain-free treatments.

Motif IR® - a safer, more efficient way of treating superficial skin texture and pore size for overall facial rejuvenation with zero downtime.



eLight™ — Safe and effective light-based elōs platform

eLight is the affordable and versatile complete facial aesthetic solution.

Customizable and upgradable - five applicators allow you to expand your aesthetic offerings at your pace.

Meet patient demand - fast and intuitive color correction and wrinkle reduction.

elure™ — Advanced skin brightening technology



- Achieve a brighter and more even skin appearance in just one month
- Developed with Melanozyme™ — a naturally occurring enzyme formulation
- Unique patented formula
- Fast, safe and gentle
- Hydroquinone free
- Clinically proven

Systems & Applications

	APPLICATIONS	PicoWay The Complete PicoSecond Platform	Alex TriVantage ¹ Q-Switched Alexandrite/ Nd:YAG	CO ₂ RE The Complete CO ₂ Solution	eLase Diode eLās Technology	eLight Light-Based eLās Technology	eLumina Advanced Bright
Body Shaping	Fat Cell Destruction						
	Non-Invasive Body Contouring						
	Cellulite Treatment						
Clinical Dermatology	Port Wine Stains		•	•			
	Bruising ²						
	Hemangiomas						
	Minor Dermal Excisions			•			
	Warts			•			
	Post Traumatic Scars ³			•			
Hair Removal	Acne Scars			•			
	Surgical Scars / Striae			•			
	Rosacea						
Pigmented Lesions	Permanent Hair Reduction				•	•	
	Beard Bumps (PFB)						
	Brown Spots	•	•			•	
	Nevus of Ota/Ito	•	•				
	Tattoo Removal	•	•				
Skin Rejuvenation ⁴	Recalcitrant Tattoos	•					
	Color Correction (Pigmented / Vascular)	•	•	•		•	
	Treat Adipose Tissue						
	Creates Elastin, Collagen & Hyaluronic Acid						
	Wrinkle Reduction			•	•	•	
	Acne Vulgaris						
	Sublative Rejuvenation™						
Vascular Lesions	triniti Plus (Total Facial Rejuvenation)						
	Light Laser Peel			•			
	Fractional Ablative Skin Resurfacing			•			
	0.2164 in			•			
	Facial and Spider Veins						
	Leg Veins						

(1) The 755nm ALEXLAZR is available in selected markets

(2) By selective targeting of injured blood vessels

(3) Via skin resurfacing

(4) Skin rejuvenation by reduction of diffuse redness and pigmentation


(5) Available as a dual wavelength 1064nm, 755nm unit or

upgradable single wavelength 1064nm and 755nm units

(6) Bipolar Fractional Radiofrequency Treatment Induces Neolastogenesis and Neocollageneesis, Hantash, Basil et al. Lasers in Surgery and Medicine, 2009

elure Advanced Skin Brightening	eIos Plus Complete Aesthetic Work Station	eTwo Sublime & Sublative RF Technology	Profound Creates Younger Skin. Profoundly.	Gentle Pro® Series Nd:YAG/ Alexandrite	Vbeam Gold Standard Pulse Dye	VelaShape III The Art of Body Contouring	UltraShape Pulsed Focused Ultrasound
							•
						•	•
						•	
					•		
				•	•		
					•		
					•		
	•	•			•		
	•	•			•		
	•			•	•		
	•			•			
	•			•	•		
				•			
	•			•	•		
			•				•
			•				
	•	•	•	•	•		
	•				•		
	•	•					
	•						
	•	•					
	•			•	•		
	•			•			

THE MOST RECOMMENDED BRAND IN THE INDUSTRY



In a Brand Perception Market Survey conducted recently by Wizer via the Harris Physician's Panel among U.S. dermatologists, Syneron Candela was named the #1 most recommended brand in the category of energy-based devices.

 SYNERON CANDELA

 Other Manufacturers



“elōs technology has been **phenomenal in my practice**. It combines energies in such a way that we can provide **exceptionally effective treatments** in the safest way possible, with the greatest level of comfort and without downtime.”

“I have owned many (Syneron) Candela devices... and have always found them to be **reliable workhorse machines** with excellent technical support.”

“Syneron Candela offers a range of **good products, good service** and have been a long-term supporter of dermatology.”

“Life long experience with (Syneron) Candela.”

INNOVATIONS THAT DEMAND ATTENTION

AWARD-WINNING TECHNOLOGY



Gentle Hair Removal

Best Laser Hair Removal:
Women's Health



UltraShape Power

Beauty Choice Award
winner: New Beauty



PicoWay

Best Picosecond Laser:
The Aesthetic Industry Awards



Profound

4 Hottest Launches in 2015:
Plastic Surgery Practice



OVER 33 BILLION MEDIA IMPRESSIONS

Breakthrough technology is bound to attract attention, which is why Syneron Candela systems have been featured in some very prominent places. Your patients are seeing these treatment options in some of their favorite magazines, on television and online.



“

Addressing slackness, crow's feet and crepey skin in a single treatment... [Profound] reaches the optimal temperature for the production of new collagen and elastin with no risk of burning.”

—Vogue, September 2015, p. 829



“

Need a tattoo gone? You'll have less pain with the PicoWay laser.”

—Women's Health, October 2015, p. 50



“

What would you say if I told you that you could get slim... without hitting the gym? UltraShape can help us do just that!”

—Good Morning America, January 5, 2016

Woman's
World

VOGUE

Modern
Aesthetics

ACCESS
HOLLYWOOD / ★

ELEVATE

Women'sHealth

Cosmetic
SURGERY TIMES

BAZAR

the
Doctors

Prevention

the
dermatologist

THE DR.
OZ
SHOW



Science. Results. Trust.



THE #1 MOST RECOMMENDED BRAND

Syneron Candela is the brand that physicians would most likely recommend to a peer.¹ That's because Syneron Candela systems are built on a solid foundation of quality and reliability. Along the way, Syneron Candela has led the way in innovations that give practices the power to deliver high-demand treatments to an ever-expanding patient base. And through it all, the Syneron Candela name has become synonymous with excellence in clinical training and technical support.

Syneron Candela is a leading global aesthetic device company with a comprehensive product portfolio and a global distribution footprint.

The Company's technology enables physicians to provide advanced solutions for a broad range of medical-aesthetic applications including body contouring, hair removal, wrinkle reduction, tattoo removal, improving the skin's appearance through the treatment of superficial benign vascular and pigmented lesions, and the treatment of acne, leg veins and cellulite.

The Company sells its products under three distinct brands, Syneron, Candela and CoolTouch, and has a wide portfolio of trusted, leading products, including UltraShape, VelaShape, GentleLase, VBeam Perfecta, PicoWay, CO₂RE, CO₂RE Intima*, Profound and elōs Plus.

Founded in 2000, the company markets, services and supports its products in 86 countries. It has offices in North America, France, Germany, Italy, Portugal, Spain, UK, Australia, China, Japan, and Hong Kong, and distributors worldwide.

SYNERON  CANDELA®

www.syneron-candela.com

¹ Proprietary Syneron Candela survey conducted by Wizer via Nielsen's online physicians panel, among dermatologists, plastic surgeons and medspas, 2015

* CO₂RE Intima is a module of CO₂RE

These products are CE approved.

©2016. All rights reserved. This marketing material may include registered trademarks, trade-names and brands in certain jurisdictions of Syneron Candela group of companies. PB91293EN